

Every Shot Counts™



Two system configurations available:





PVX100 Dual Configuration

Simultaneous administration of **two vaccines**.





PVX100 System Overview





A comprehensive solution



Mobile App

Seamless programming of multiple systems with pre-defined vaccination plans.

Real-time process oversight and injection quality control with alerts or vaccination issues (including missed injections, needle replacement, incorrect vaccine temperature, and more.

Offline capabilities for rural areas with no internet connectivity.

Web App

Stores vaccination data in a secure, password-protected
environment, remotely accessible 24/7 from the cloud.

- Analyze and compate historical data records to implement data-driven decisions for continuous improvement.
- Advanced Business Intelligence dashboards and reports for intuitive data interpretation.







System Specifications

System Configuration	Single & Dual				Dual	
Needles	1 needle				2 needles	
Vaccine sources	1-2 vaccines				2 vaccines	
Route of administration	IM & SC				IM	
Injection Head	Regular	Slim	Slim Long*	Slim XL*	Twin	Twin Slim*
Needle length in inches	³ / ₁₆ , ¹ / ₄ , ³ / ₈ , ¹ / ₂	3⁄8	³ /8, ¹ /2	⁵ /8, ³ /4	¹ ⁄4, ³ ⁄8, ¹ ⁄2	1⁄2, ⁵ ⁄8
Needle length in mllimeters	5, 6.5, 9.5, 12.7	9.5	9.5, 12.7	15.9, 19	6.5, 9.5, 12.7	12.7, 15.9
Max. needle penetration	3-9mm	3-6mm	3-9mm	9-15.5mm	3-9mm	3-9mm
Needle hub	plastic / metal	DIASTIC ONIV				
Needle gauge Ø	18G, 19G, 20G (optional thinner needle with low dosage)					
Dosage range in milliliters	0.2 - 1.0ml (higher/lower dosage optional)					
Emulsion compatibility	oil & water-based vaccines					

* available soon



For inquiries on the pHi-Tech PVX100 system, please reach out to your local Phibro representative.

For additional information visit our website **phitech.pahc.com**

©2024 pHi-Tech Animal Health Technologies LTD

Every Shot Counts[™] is a trademark owned by or licensed to Phibro Animal Health Corporation or its affiliates.